STRATTICE™ Reconstructive Tissue Matrix  
Perforated

Instructions  
For Use

Surgical Mesh

# Table of Contents

## Sommaire

Instructions for Use   English.....	3–6
Instructions d'utilisation   Français.....	7–11

## Device Description

STRATTICE™ Reconstructive Tissue Matrix Perforated (STRATTICE™ TM or STRATTICE™ surgical mesh) is a surgical mesh that is derived from porcine skin and is processed and preserved in a patented phosphate buffered aqueous solution containing matrix stabilizers. STRATTICE™ Tissue Matrix is intended to perform as a surgical mesh for soft tissue repair while presenting a scaffold to the patient's tissue. The structural properties minimize tissue adhesions to the mesh. The surgical mesh consists of a sterilized sheet of processed porcine dermis provided in prescribed geometric configurations and packaged in a double pouch.

Use of the surgical mesh provides for a strong and biocompatible implant and will incorporate into the patient's tissue with associated cell and microvascular ingrowth.

Animal studies show a low incidence of adhesion to the STRATTICE™ Tissue Matrix based on observation of minimal visceral tissue attachment.

## Indications

STRATTICE™ Reconstructive Tissue Matrix Perforated is intended for use as a soft tissue patch to reinforce soft tissue where weakness exists and for the surgical repair of damaged or ruptured soft tissue membranes. Indications for use include the repair of hernias and/or body wall defects which require the use of reinforcing or bridging material to obtain the desired surgical outcome.

STRATTICE™ Reconstructive Tissue Matrix Perforated is supplied sterile and is intended for single patient one-time use.

## Contraindications

- This surgical mesh is derived from a porcine source and should not be used in patients with known sensitivity to porcine material.
- Polysorbate 20 is a component of the aqueous phosphate buffered solution and therefore STRATTICE™ Tissue Matrix should not be used in patients with a known sensitivity to this material.

## Warnings

- **Do not resterilize.**
- Do not use if the package is opened or damaged. Do not use if seal is broken or compromised.
- Do not use if the temperature monitoring device does not display "OK".
- After use, handle and dispose of all unused product and packaging in accordance with accepted medical practices and applicable national and regional environmental laws on disposal of packaging and biological materials.
- STRATTICE™ Tissue Matrix cannot be reused once it has been removed from the packaging and/or is in contact with a patient without increased risk of patient-to-patient contamination and subsequent infection.

## Precautions

- Discard surgical mesh if handling has caused possible damage or contamination.
- Discard surgical mesh if it is past the use-by-date of the product (indicated as 4 digit year, 2 digit month, and 2 digit day [YYYY-MM-DD]).
- Ensure that the surgical mesh is put into a sterile basin and covered with room temperature sterile saline or sterile Ringer's lactate for a minimum of 2 minutes prior to implantation.

# STRATTICE™

## Instructions For Use | English

- Place surgical mesh in maximum possible contact with healthy, well-vascularized tissue to promote cell ingrowth and tissue remodeling.
- The surgical mesh should be hydrated and moist when the package is opened. If the surgical mesh is dry, do not use.
- In significantly contaminated or infected cases, utilize bioburden-reducing techniques to minimize contamination levels at the surgical site, including, but not limited to, appropriate drainage, debridement, negative pressure therapy, and/or antimicrobial therapy prior and in addition to implantation of the surgical mesh. (Presence of a significant microbial load may impact overall performance of the surgical mesh.)
- Large hernia defects and a bridging mesh technique are risk factors for hernia recurrence. Likewise, in large abdominal wall defect cases where midline fascial closure cannot be obtained, with or without separation of components techniques, utilization of STRATTICE™ TM in a bridged fashion is associated with a higher risk of hernia recurrence than when used to reinforce fascial closure.
- If a tissue punch-out piece is visible, remove using aseptic technique before implantation.
- Certain considerations should be made in order to reduce the risk of adverse events when performing surgical procedures using STRATTICE™ TM. Please see the following sections for more information: Product/Patient Selection, Technique Guidance, and Post-Operative Care.

### Adverse Events

Potential adverse events are those typically associated with surgical mesh materials and their implantation procedures including, but not limited to, infection, foreign body response, hematoma, seroma formation, failure to integrate, recurrence of tissue defect, bulging, fistula formation, tissue necrosis, inflammation and adhesion formation. If an unanticipated event occurs, alteration of surgical plan may be necessary at the surgeon's discretion.

### Storage

- The surgical mesh is a sterile medical device that should be stored in a clean, dry location at -8°C to 30°C.
- Refer to the temperature monitor located on the product carton to ensure that product has been stored within its temperature limits. Only use the product if the included temperature monitor displays "OK" on the screen. If screen displays anything other than "OK", do not use the product.
- It is to be stored in its original packaging.
- The use-by-date of the product is indicated as 4 digit year, 2 digit month, and 2 digit day (YYYY-MM-DD).

### Sterilization

This product is sterilized by E-Beam.

### Product/Patient Selection

The following considerations are intended to serve only as general guidelines. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care. Physicians should be familiar with surgical procedures and techniques involving biological surgical mesh before using STRATTICE™ TM.

- Carefully consider the risk/benefit balance of using STRATTICE™ TM in patients with significant comorbidities; including but not limited to, obesity, smoking, diabetes, immunosuppression, malnourishment, poor tissue oxygenation (such as COPD), and pre- or post-operative radiation.

## Instructions For Preparing STRATTICE™ Tissue Matrix For Surgical Use

These instructions are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care. Physicians should be familiar with surgical procedures and techniques involving biological surgical mesh before using STRATTICE™ Tissue Matrix.

### Required Materials

- Sterile forceps
- Soaking fluid: room temperature sterile saline or room temperature sterile lactated Ringer's solution.
- One sterile basin per piece of surgical mesh

### Preparation

1. Open the carton and remove the pouch.
2. Peel open the outer pouch and remove the inner pouch using aseptic technique. The inner pouch is sterile and may be placed directly into the sterile field.
3. Open the inner pouch carefully and aseptically remove the surgical mesh. Always use sterile gloves or forceps when handling the surgical mesh.
4. Place the surgical mesh in a sterile basin and submerge in room temperature sterile saline or room temperature sterile lactated Ringer's solution for a minimum of 2 minutes.
5. Keep the surgical mesh in the room temperature sterile solution until ready for implantation. Surgical mesh can be retained in sterile solution for a maximum of 4 hours.

### Implantation

1. Prepare the surgical site using standard techniques.
2. STRATTICE™ Tissue Matrix may be folded, trimmed or cut as required to fit the surgical site using aseptic technique, ensuring allowance for overlap.
3. Transfer the surgical mesh to the surgical site using sterile gloves or forceps.
4. Suture the surgical mesh into place.
5. Complete the standard surgical procedure.
6. Discard any unused portions of the surgical mesh as per institutional procedures.

The following suggestions are not intended to supersede institutional protocols or professional clinical judgment concerning patient care. Physicians should be familiar with surgical procedures and techniques involving biological surgical mesh before using STRATTICE™ TM.

### Technique Guidance

- As with any surgical implant, careful aseptic technique should be practiced and contact of the mesh with the patient's skin should be minimized.

# STRATTICE™


Instructions For Use | English

- Tension and suture placement are application dependent. For hernia repair applications, surgical experience with biological surgical mesh indicates that suturing the STRATTICE™ TM under physiologic tension with a minimum of 3cm–5cm overlap, or as much as required to reach healthy adjacent tissues, may produce improved outcomes.
- Use of permanent sutures is recommended.

## Post-Operative Care

- Proper post-operative care may include the use of appropriate drainage and negative pressure therapy.
- As with any post-operative care, aseptic technique should be practiced when required to minimize contamination to the surgical wound.

## Definitions

	PERFORATED		OPEN HERE		CONTENTS
	BATCH CODE		INNER POUCH		OUTER POUCH
	DO NOT USE IF PACKAGE IS DAMAGED		CORRECT SETTING		INCORRECT SETTING
	FOR SINGLE USE ONLY		DO NOT RESTERILIZE		ATTENTION SEE INSTRUCTIONS FOR USE
	CATALOG NUMBER		USE BY DATE		TEMPERATURE LIMITATIONS
	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY		STERILIZED USING IRRADIATION		MANUFACTURED BY

**Rx ONLY** CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

## Customer Support

For product complaints and potential adverse events, please contact your local Sales Representative, or 1-800-433-8871.

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Allergan® and its design are trademarks of Allergan, Inc.
STRATTICE™ is a trademark of LifeCell Corporation, an Allergan affiliate.
© 2018 Allergan. All rights reserved.
LifeCell Corporation